


Zapp is certified to ISO 9001



Grade Zapp® 1.4108

Zapp® 1.4108 is a nitrogen-alloyed, high corrosion resistant martensitic steel with excessive toughness at hardness up to 60 HRC. Partly replacing carbon with nitrogen results in a much higher corrosion resistance and wear resistance compared to conventional, hardenable martensitic grades. In combining the Pressure Electro-Slag-Remelting-Process (PESR) with an elaborate forging technique, an extremely high purity level of a fine and homogeneous microstructure can be achieved. This implies excellent machinability, outstanding polishing and high dimensional stability after heat treatment. Consequently, Zapp® 1.4108 is the ideal grade for bending-stressed or break-endangered medical instruments which are in contact with highly corrosive mediums.

Typical Fields of Applications

- Medical Instruments e. g.
- Drills
- Screwdrivers
- Chisels
- Saw Blades, Cutting Tools

Corresponding Standards

DIN X30CrMoN15-1



UNS S42027 acc. ASTM F899

Microstructure Pictures

In comparison to conventional martensitic grades, Zapp® 1.4108 shows significantly smaller primary carbides.

Zapp® 1.4108

Zapp® 1.4112 (AISI 440 B)



Typical Chemical Composition*

C	Si	Mn	Cr	Mo	N
0.30	0.60	0.40	15.00	1.00	0.40

* average in mass-%

Product Conditions*

Bars, drawn, straightened, ground	Tensile [MPa]	700-900
-----------------------------------	---------------	---------

* Special conditions on request

Physical Properties

Modulus of Elasticity E 20 °C [GPa]	223
Specific Gravity [kg/dm³]	7.72
Thermal Conductivity 20°C [W/m*K]	14
Coefficient of Thermal Expansion α [10 ⁻⁶ *K ⁻¹]	
20 – 100 °C	10.4
20 – 200 °C	10.8
20 – 300 °C	11.2
20 – 400 °C	11.6
20 – 500 °C	11.9
Specific Heat c at 20 °C [kJ/kg*°C]	430
Electric Resistivity ρ at 20 °C [Ω*mm²/m]	0.8

Heat Treatment

Soft Annealing

780 – 820 °C/ 7 h / Cooling: furnace or air

Stress Relief Annealing

150 – 220 °C/ 2 x 2 h/ Cooling: Air

Hardening

1,000 – 1,050 °C/ 0,5 h/ Cooling: Oil

Hardening has to be conducted under nitrogen partial pressure to prevent reduction or increase of the nitrogen content.


Tempering

100 – 475 °C/ 2 x 2 h/ Cooling: Air

Subzero Refrigeration

-80 – -196 °C/ 1 h/ applied to eliminate remaining austenite at hardening temperatures of > 1,010 °C

Tempering Chart (Hardening with Subzero Refrigeration)



Corrosion Resistance

Through the addition of nitrogen, Zapp® 1.4108 shows an exceptional corrosion resistance.

Hot Working

Forging at 1,220 – 1,000 °C

Magnetism

Zapp® 1.4108 is magnetizable.

Machining


Zapp® 1.4108 is characterized by an outstanding machinability.

Polishability

Zapp® 1.4108 shows excellent abilities for grinding and polishing.

Comparison of Corrosion Resistance

PREN (Pitting Resistance Equivalence Number)



Corrosion Resistance of Zapp® 1.4108 in comparison to conventionally hardenable martensitic grades.

MEDICAL ALLOYS

Zapp Precision Metals GmbH

Letmather Straße 69

58239 Schwerte

medicalalloys@zapp.com

www.zapp.com

Further information regarding our products and locations are available in our image brochure and under www.zapp.com

The illustrations, drawings, dimensional and weight data and other information included in this data sheet are intended only for the purposes of describing our products and represent non-binding average values. They do not constitute quality data, nor can they be used as the basis for any guarantee of quality or durability. The applications presented serve only as illustrations and can be construed neither as quality data nor as a guarantee in relation to the suitability of the material. This cannot substitute for comprehensive consultation on the selection of our products and on their use in a specific application. The brochure is not subject to change control.

Last revision: June 2026